

# THE ASCLEPION EFFECT

## QuadroStar+ Lipo

High power laser  
assisted liposuction

faster, safer & flexible

FAT REMOVAL

SLIMMING

BODY CONTOURING

SKIN TIGHTENING

COLLAGEN REMODELLING



Technology made in Germany

**Asclepi<sup>o</sup>n**  
Laser Technologies

# QuadroStar<sup>+</sup> Lipo

## USPs/Benefits:

### QuadroStar<sup>+</sup> 980 technical advantages

- proven technology
- table - top system, transportable
- best price/value ration
- complete Lipolysis-Sets with multiple fiber sizes
- very fine, easy & quick exchangeable cannulas
- luer lock adapter - easy & safe
- automatic fiber size detection - safe & easy handling
- 532 upgrade possible




Before treatment



1 month after treatment

### QuadroStar<sup>+</sup> 980 application advantages

- QuadroStar 980 laser assisted lipolysis - multiple effects
- volume reduction and skin tightening
- smoother skin with less deformities compared to conventional liposuction
- faster healing with smaller incisions, less bruising, reduced pain minimal blood loss and swelling
- less down - time
- minimal traumatization due to very fine cannulas
- easy application on critical locations - more flexible than liposuction



Handpiece with 5 cannulas

### QuadroStar<sup>+</sup> 980 specifications

The first table top, high power cw - mode laser for lipolysis

- 25 Watt on tissue
- 980 nm
- cw, varix, repeat and single mode
- upgradeable with second wavelength (532 nm)
- 19 kg - portable

### Lipo Set

- 2 handpieces
- 2 fibers (300  $\mu$ m and 600  $\mu$ m)
  - 3 cannulas for 300  $\mu$ m, length 50, 100, 150 mm
  - 2 cannulas for 600  $\mu$ m, length 80, 150 mm